

5 ideas to inspire communication using Message Window

Message Window lets you build complete phrases on a Lingraphica device. Instead of selecting one card at a time in Talk, you can combine words into full sentences — then tap “Speak” to say it all at once. This makes it easier to share needs, stories, and conversations.



There's not one "right" way to use Message Window. Whether you're a device user, a loved one, or a speech-language pathologist — here are a few quick tips to keep in mind as you explore this feature:



Start with words you use every day. Create messages that relate to daily routines, interests, and needs. For example, you might say “I’d like a triple espresso” or “I only drink decaf tea” — depending on your unique preferences.



Personalize your AAC device. Include images and words that are special to you — like a picture of your pet and their name or your favorite dish at a local restaurant. This makes it easier to use your device in real life (and share about the things you love the most).




Get free support whenever you need it. Remember, you're never on your own. Lingraphica offers free, personalized training for device users, families, and SLPs. We'll help you get comfortable with Message Window and every part of your device. Learn more at www.lingraphica.com/support.

Your Lingraphica device will help you go beyond single words, so you can say exactly what's on your mind. Next, you'll see 5 different ideas that you can use as inspiration when you communicate with Message Window.

Note: Some messages you see below may come pre-programmed on the device, while others will need to be added. Our team can help you with adding personalized new messages to your Lingraphica device during or after your free device trial.



1. Start with an introduction.



This is a card.

This is a folder.

It has many cards inside it. Folders help you share even more details about your favorite hobbies.




Tip:
Add common questions to keep the conversation going.

2. Share how you feel.



This is a "web card."
When you tap this one, it will take you to the "Wheel of Emotions" website. You can create web cards that link to any website.




Communicating about pain
can be difficult, but this folder will make it simpler. It includes a pain scale, numbers 0-10, and cards for each body part.

3. Talk about your favorite things.



Inside this folder,
this Lingraphica user has
pictures of their favorite players,
links to their favorite football
websites, and more. This helps
them talk about their favorite
team in more depth.




**You can add any picture
from the internet**
to match your message.

4. Advocate for yourself and your needs.



“Carrier phrases”
like “I need to” can help
you start many different
statements.




Each of these messages
could flow easily after the
carrier phrase, “I need to.”

5. Play games and have fun.



This card features a “dynamic image.”

When you select this card, you’ll see a short video of the boy catching the ball. Dynamic images can make communication more relatable and intuitive for people with language challenges.



AAC isn’t just for expressing wants and needs.

Include phrases to comment on a variety of situations.



Remember:

These phrases are just a few ideas to get started! With Message Window, you can build complete thoughts and say exactly what’s on your mind — one phrase at a time.

Need help getting started with a Lingraphica device?
We’re here for you.

Contact us at 866-901-4886 or lingraphica.com/contact-us/.